



Meet the new face of performance

Tissue Cassettes from Epredia



The entire range of cassettes from EpreDia™ has been redesigned for better performance – from printing through laboratory processes to archiving

Printing

- The all-new surface finish and reinforced print faces are engineered to optimize print quality when used with PrintMate™ AS cassette printers from EpreDia
- All-new PrintMate AS-compatible tubes are precisely engineered to reduce jamming and prevent downtime

Grossing, Embedding and Processing

- Opening and closing forces have been carefully adjusted for ease of use, while safeguarding valuable specimens during processing
- Design and flow characteristics tested to ensure efficient processing

Sectioning

- The new cassettes are precision engineered for a solid connection to EpreDia microtomes





The mark of performance

It's easy to identify cassettes designed for performance. Just look for the unique PrintMate Certified icon on the tab.





**PrintMate Slotted
Cassette from Epredia**

- Pore style: Slotted
- For routine specimens
- Lid hinge: Reversible

B851090 □□ 1,000 cassettes,
loose, with loose-packed lids

B851120 □□ 750 cassettes
10 tubes of 75 cassettes,
with loose-packed lids



**Thinline
Slotted Cassette**

- Pore style: Slotted
- For routine specimens
- Lid hinge: Front edge

B851736 □□ 750 cassettes
10 tubes of 75 cassettes, with
color-matched, attached lids



**PrintMate Biopsy
Cassette from Epredia**

- Pore style: Small squares
- For biopsy specimens
- Lid hinge: Reversible

B851080 □□ 1,000 cassettes,
loose, with loose-packed lids

B851110 □□ 750 cassettes
10 tubes of 75 cassettes,
with loose-packed lids



**Thinline
Biopsy Cassette**

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Front edge

B851737 □□ 750 cassettes
10 tubes of 75 cassettes, with
color-matched, attached lids



**PrintMate Embedding
Cassette from Epredia**

- Pore style: Round holes
- For routine specimens
- Lid hinge: Reversible

B851040 □□ 1,000 cassettes, loose,
with no lids – for use with metal lids

B851010 □□ 750 cassettes, 10 tubes of 75
cassettes, with no lids – for use with metal lids

B851070 □□ 1,000 cassettes,
loose, with loose-packed lids

B851100 □□ 750 cassettes, 10 tubes
of 75 cassettes, with loose-packed lids

1001003 Metal lids for cassettes





Slotted Cassette 2

- Pore style: Slotted
- For routine specimens
- Lid hinge: Front edge

B851739 □□ 1,000 cassettes, loose, with loose-packed lids

B851729 □□ 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



TissueLoc™ Cassette (Cassette IV) from EpreDia

- Pore style: Slotted
- For routine specimens
- Lid hinge: Back edge
- Large tab for easy opening

B851954 □□ 1,000 cassettes, loose, with loose-packed lids

B851720 □□ 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



HistoScreen™ Cassette from EpreDia

- Pore style: Mesh
- For very small biopsy specimens
- Eliminates need for sponges or bags
- Lid hinge: Back edge

B851000 □□ 1,000 cassettes, loose, with loose-packed lids

B851001 □□ 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



Biopsy Cassette 2

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Front edge

B851738 □□ 1,000 cassettes, loose, with loose-packed lids

B851730 □□ 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



TissueLoc™ Biopsy Cassette (Cassette IV) from EpreDia

- Pore style: Small square
- For biopsy specimens
- Lid hinge: Back edge
- Large tab for easy opening

B851925 □□ 1,000 cassettes, loose, with loose-packed lids

B851750 □□ 750 cassettes, 10 tubes of 75 cassettes, with loose-packed lids



Fine Pore Cassette

- Pore style: Small holes
- For very small biopsy specimens
- Lid hinge: Front edge

B851735 □□ 1,000 cassettes, loose, with color-matched, attached lids

B851734 □□ 750 cassettes, 10 tubes of 75 cassettes, with color-matched, attached lids

Ordering is easy.

Simply select your cassette and configuration number above, add the two digit color code from below.

Example: Light Green TissueLoc Biopsy Cassette in Tubes = B851925LG



WH
White



GY
Grey



TN
Tan



BE
Beige



YW
Yellow



GN
Green



LG
Light Green



LA
Lavender



BL
Blue



LI
Lilac



PK
Pink



OR
Orange



PrintMate AS on-demand cassette printer

Our cassettes are precisely engineered for optimal performance when paired with the PrintMate™ AS cassette printer from EpreDia.

Eliminate waste and streamline your labeling process by printing cassettes on-demand with the PrintMate AS at the grossing station.

- Save time – prints each cassette in ten seconds or less
- Print only the cassettes you need, right at the point of use
- Print barcodes on cassettes to automate downstream processes
- Eliminate errors endemic to handwritten processes
- Suitable for both on-demand and batch printing
- Optional collection system for efficient and organized batch printing
- Free your workspace, the PrintMate printer's compact footprint fits easily and integrates into your laboratory's workflow



Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2023 EpreDia Group. All rights reserved. All trademarks are the property of EpreDia Group and its subsidiaries unless otherwise specified. **M44047 R0523**

