


Clear, borosilicate glass for coverslipping slides in the lab

Epredia Coverglass

Available in rectangles, squares, circles, Superslip™ and premium options





Epredia Coverglass

- Coverslips are manufactured with colorless borosilicate glass. Epredia coverslips feature high spectral-transmission quality and a finely adapted refractive index suited for microscopy
- They are available in a variety of sizes, thicknesses and shapes
- Thicknesses range from
 - #0 (0.085 to 0.13 mm)
 - #1 (0.13 to 0.16 mm)
 - #1.5 (0.16 to 0.19 mm)
 - #2 (0.19 to 0.23)
- All coverglass is sold by the case which includes 10 1oz boxes of coverglass.

Superslip Coverglass


Superslip coverslips provide a more slippery surface than traditional coverglass to enhance usability and compatibility with automated coverslippers.

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
24X40-1W-YO-001G	24 x 40 mm	#1	82
24X50-1W-YO-001G	24 x 50 mm	#1	54
152450S	24 x 50 mm	#1.5	48
12455S	24 x 55 mm	#1	48
24X60-1W-YO-001G	24 x 60 mm	#1	65

Premium Coverglass

Premium coverglass from Epredia is packaged in specialty boxes and includes a desiccant

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
22X40-1-002G	22 x 40 mm	#1	89
22X22-1-002G	22 x 22 mm	#1	162
24X40-1-002G	24 x 40 mm	#1	82
24X50-1-002G	24 x 50 mm	#1	65
24X60-1-002G	24 x 60 mm	#1	54
48X60-1-002G	48 x 60 mm	#1	27



Coverglass Rectangles

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
22X30-1-001G	22 x 30 mm	#1	119
22X30-1.5-001G	22 x 30 mm	#1.5	103
22X32-1.5-001G	22 x 32 mm	#1.5	90
22X40-1-001G	22 x 40 mm	#1	89
22X40-1.5-001G	22 x 40 mm	#1.5	72
22X40-2-001G	22 x 40 mm	#2	59
22X50-1-001G	22 x 50 mm	#1	71
22X50-1.5-001G	22 x 50 mm	#1.5	58
22X60-1-001G	22 x 60 mm	#1	58
22X60-1.5-001G	22 x 60 mm	#1.5	48
24X30-1-001G	24 x 30 mm	#1	109
24X30-1.5-001G	24 x 30 mm	#1.5	94
24X40-1-001G	24 x 40 mm	#1	82
24X40-1.5-001G	24 x 40 mm	#1.5	66
24X50-1-001G	24 x 50 mm	#1	65
24X50-1.5-001G	24 x 50 mm	#1.5	56
24X50-2-001G	24 x 50 mm	#2	43
24X55-1-001G	24 x 55 mm	#1	59
152455	24 x 55 mm	#1.5	51
24X60-1-001G	24 x 60 mm	#1	54
24X60-1.5-001G	24 x 60 mm	#1.5	47
24X60-2-001G	24 x 60 mm	#2	36
25X60I-2-4789-1-LS Lifter Slip	25 x 60 mm	N/A	39



1oz package of coverglass




Package of 10 1oz boxes of coverglass

Coverglass Squares

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
18X18-1-001G	18 x 18 mm	#1	242
18X18-1.5-001G	18 x 18 mm	#1.5	194
22X22-0-001G	22 x 22 mm	#0	240
22X22-1-001G	22 x 22 mm	#1	162
22X22-1.5-001G	22 x 22 mm	#1.5	130
22X22-2-001G	22 x 22 mm	#2	120
24X24-1-001G	24 x 24 mm	#1	136
25X25-1-001G	25 x 25 mm	#1	121
25X25-2-001G	25 x 25 mm	#2	92

Coverglass Circles

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
10CIR-1-001G	10 mm	#1	962
12CIR-1-001G	12 mm	#1	685
12CIR-1.5-001G	12 mm	#1.5	567
12CIR-2-001G	12 mm	#2	467
15CIR-1-001G	15 mm	#1	460
15CIR-1.5-001G	15 mm	#1.5	362
15CIR-2-001G	15 mm	#2	307
18CIR-1-001G	18 mm	#1	320
18CIR-2-001G	18 mm	#2	213
18CIR-1.5-001G	18 mm	#1.5	252
22CIR-1-001G	22 mm	#1	199
22CIR-2-001G	22 mm	#2	143
22CIR-1.5-001G	22 mm	#1.5	168
25CIR-1-001G	25 mm	#1	165
25CIR-1.5-001G	25 mm	#1.5	130
25CIR-2-001G	25 mm	#2	110
30CIR-1.5-001G	30 mm	#1.5	91
40CIR-1-001G	40 mm	#1	60
50CIR-1.5-001G	50 mm	#1.5	33



Coverglass Hoppers for ClearVue™ Automated Coverslipper

- To be used with Epreidia ClearVue™ Automated Coverslipper
- 2 hoppers per package

Cat. Number	Dimensions	Thickness	Approx number of coverslips per package
A79210166	24x50mm	#1	575
A79210050	24x40mm	#1.5	500
A79210051	24x50mm,	#1.5	500
A79210052	24x55mm	#1.5	500

To learn more about Epreidia microscope slides, coverglass and automated coverslipping solutions, visit [epredia.com](https://www.epredia.com)

