



Epredia Microscope Slides and Coverglass Brochure




	Histology	Cytology	Microbiology	Hematology	Urinalysis	Other
Superfrost PLUS™	✓	✓				
Superfrost PLUS Control	✓	✓				
Superfrost PLUS Gold	✓	✓				
Colorfrost PLUS™	✓	✓				
SlideMate PLUS™	✓	✓				
SlideMate™ Laser PLUS	✓	✓				
Polysine™ Adhesion	✓	✓				
Permafrost™ Adhesion	✓	✓				
Inkjet PLUS	✓	✓				
Colormark PLUS™	✓	✓				
Superfrost™		✓	✓	✓	✓	
Superfrost™ Control		✓	✓	✓	✓	
Colorfrost™		✓	✓	✓	✓	
SlideMate™		✓	✓	✓	✓	
SlideMate™ Laser		✓	✓	✓	✓	
Permafrost™		✓	✓	✓	✓	
Colormark™		✓	✓	✓	✓	
Premium Plain and Economy Plain			✓	✓	✓	
Premium Frosted and Frosted			✓	✓	✓	
Specialty						✓
Cytology		✓				
Petrographic						✓
Microarray						✓

All slides feature

- Dimensions 25 x 75 x 1.0 mm unless otherwise stated
- 90°, finely ground edges unless otherwise stated
- 3/4" tab unless otherwise stated

Microscope slides may be sold by the gross, 144 slides, 1/2 gross, 72 slides, or other configurations. Slide quantities and purchasing information can be found in the header of each table.



Microscope Slides

Epredia microscope slides are optically clear, exceptionally flat, and made only from low fluorescence glass. High quality glass is critical in many sections of the laboratory.

Based on the properties of each slide type, they are more often used in a specific area of the laboratory, as shown in the chart on the left.



Adhesion Slides

- Slides secure frozen and formalin-fixed paraffin embedded tissue sections as well as cytology preparations
- Tissue loss during staining can be minimized thereby resulting in rapid turn-around times and reduced risk associated with tissue loss and repeat testing




Non-Adhesion Slides

- Manufactured from extra-white soda-lime glass with very low iron content
- Clean, flat and ready to use immediately
- 25 x 75 x 1.0 mm unless otherwise stated
- 1mm thickness unless otherwise stated
- Slides labeled for export are packaged in a moisture control barrier bag

- Clipped corner slides feature a 45° edge
- Clipped corner slides use the same adhesion technology as standard 90° corner slides for each family



Standard corners



Clipped corners

Superfrost PLUS Slides

- Superfrost PLUS slides are available with a white tab
- Superfrost PLUS and Colorfrost PLUS adhesion slides are designed to secure tissue sections without the need for laboratory-applied coatings

Epredia Superfrost PLUS Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Superfrost PLUS Slides	4951PLUS-001
Superfrost PLUS Slides with EPREDIA imprint	6776214
w/ clipped corners	4951PLUS-4-001
w/ 1" tab	4941PLUS-2681-006

- Slides packaged in a moisture-control barrier bag

Epredia Superfrost PLUS Slides
 - Unit of Measure - Case (CV)
 - 2,880 slides, packaged as 2 cases, each with 20 boxes of 72 slides

Description	POR Cat. Number
Superfrost PLUS Slides (for export)	4951PLUS-001E
w/ clipped corners (for export)	4951PLUS-4-001E

Superfrost slides

- Easy to write and print on while also readable and durable
- May be used in general laboratory applications including histology, cytology, hematology, microbiology and urinalysis
- Printed tab is impervious to virtually all common laboratory chemicals and reagents




Epredia Superfrost Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Superfrost Slides	4951-001T
w/ Clipped Corner	4951-4-001
3 x 1"	4951F-001T
Cygen™ Slides	4951-7140-006
Cytogenetics Slides	CGN-4951-006

- Slides packaged in a moisture-control barrier bag

Epredia Superfrost Slides
 - Unit of Measure - Case (CV)
 - 2,880 slides, packaged as 2 cases, each with 20 boxes of 72 slides

Description	POR Cat. Number
Superfrost Slides (for export)	4951-001ET



Colorfrost PLUS slides

- Same adhesion properties as Superfrost PLUS
- Colorfrost PLUS slides are available in a variety of colors

Epredia Colorfrost PLUS Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Blue	5951PLUS-001
Green	6951PLUS-001
Pink	7951PLUS-001
Yellow	8951PLUS-001
Aqua	9951APLUS-001
Burnt Orange	9951BOPLUS-001
Gold	9951GLPLUS-001
Lilac	9951LPLUS-001
Tan	9951TPLUS-001

Colorfrost Slides

- Same technology as Superfrost slides
- Available in a variety of colors

Epredia Colorfrost Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Blue	5951-001T
Green	6951-001T
Pink	7951-001T
Yellow	8951-001T
Lilac	9951L-001

- Slides packaged in a moisture-control barrier bag

Epredia Colorfrost Slides
 - Unit of Measure - Case (CV)
 - 2,880 slides, packaged as 2 cases, each with 20 boxes of 72 slides

Description	POR Cat. Number
Blue, 20 gross (for export)	5951-001ET
Green, 20 gross (for export)	6951-001ET
Pink, 20 gross (for export)	7951-001ET
Yellow, 20 gross (for export)	8951-001ET




Colorfrost PLUS slides

- Same adhesion properties as Superfrost PLUS
- Colorfrost PLUS slides are available in a variety of colors

Epredia Colorfrost PLUS Slides with EPREDIA imprint
 - Unit of Measure - Case (CAR)
 - 1,000 slides, packaged as 10 boxes of 100 slides

Description	KZOO Cat. Number
Blue	9991001
Green	9991003
Pink	9991004
Yellow	9991002
Aqua	9991013
Burnt Orange	9991011
Gold	9991014
Lilac	9991009
Tan	9991012



SlideMate PLUS and SlideMate Slides

- Formulated for use with SlideMate AS and SlideMate Pro printers
- Self-contained packages ready for loading into the SlideMate hoppers
- Simply remove the cover from the slide package and drop the sleeve into the hopper
- Available with PLUS formula to provide the maximum tissue adhesion

Eprelia SlideMate PLUS Slides


- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR and KZOO Cat. Number
White	TT-4041I-PS-8218-1
Blue	TT-5041I-PS-8218-1
Green	TT-6041I-PS-8218-1
Pink	TT-7041I-PS-8218-1
Yellow	TT-8041I-PS-8218-1

Eprelia SlideMate Slides

- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR and KZOO Cat. Number
White	TT-4041I-8217-001
Blue	TT-5041I-8217-001
Green	TT-6041I-8217-001
Pink	TT-7041I-8217-001
Yellow	TT-8041I-8217-001



SlideMate Laser PLUS and SlideMate Laser Slides

- Formulated to be used with SlideMate Laser printer
- Self-contained packages ready for loading into the SlideMate hoppers
- Simply remove the cover from the slide package and drop the sleeve into the hopper
- Available with PLUS formula to provide the maximum tissue adhesion

Eprelia SlideMate Laser PLUS Slides


- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR and KZOO Cat. Number
White	LS-4041I-PS8523-1
Blue	LS-5041I-PS8523-1
Green	LS-6041I-PS8523-1
Pink	LS-7041I-PS8523-1
Yellow	LS-8041I-PS8523-1

Eprelia SlideMate Slides

- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR and KZOO Cat. Number
White	LS-4041I-8524-1
Blue	LS-5041I-8524-1
Green	LS-6041I-8524-1
Pink	LS-7041I-8524-1
Yellow	LS-8041I-8524-1



SlideMate AS
Cat. Number A83900005



SlideMate Pro
Cat. Number A83900015




SlideMate Laser
Cat. Number 8000048A

Permafrost Adhesion Slides

- Secures frozen and formalin-fixed paraffin embedded tissue sections as well as cytology preparations


Epredia Permafrost Adhesion Slides		
- Unit of Measure - Case (CV)		
- 1,440 slides, packaged as 20 boxes of 72 slides		
Description	POR Cat. Number	KZOO Cat. Number
White	640-ADH-001	640-ADH
Blue	642-ADH-001	N/A
Green	643-ADH-001	N/A
Pink	644-ADH-001	N/A
Yellow	646-ADH-001	N/A
Burnt Orange	940BO-ADH-001	N/A
Tan	940T-ADH-001	N/A



Permafrost Slides

- General laboratory applications slides
- Frosted end for marking slides

Epredia Permafrost Slides		
- Unit of Measure - Case (CV)		
- 1,440 slides, packaged as 20 boxes of 72 slides		
Description	POR Cat. Number	KZOO Cat. Number
White	640-001T	4445
White, Clipped corners	640-4-7431-001	N/A
Aqua	N/A	4441
Blue	5951-001T	4446
Gold	N/A	4442
Green	6951-001T	4450
Lilac	9951L-001	4444
Burnt Orange	N/A	4448
Pink	7951-001T	4451
Tan	N/A	4440
Yellow	8951-001T	4447
Red	N/A	4443




Polysine Slides

- Adhesion coating derived of the chemistry for poly-L-Lysine attracts paraffin-embedded sections to slide surface

Epredia Polysine Slides	
- Unit of Measure - Case (CV)	
- 1,440 slides, packaged as 20 boxes of 72 slides	
Description	POR Cat. Number
White Tab	P4981-001

Epredia Polysine Slides	
- Unit of Measure - Case (CAR)	
- 1,440 slides, packaged as 20 boxes of 72 slides	
Description	KZOO Cat. Number
White Tab	6776215
Blue Tab	6776216



Superfrost PLUS Gold adhesion slides

- Complex adhesion technology which attracts and chemically bonds fresh or formalin-fixed frozen tissue sections firmly to the surface of the slide
- Recommended for hard-to-hold tissue samples like bone, brain, breast and skin
- Ideal for special stains, immunocytochemical and in-situ DNA hybridization on these tissue types

Epredia Superfrost PLUS Gold Slides	
- Unit of Measure - Case (CV)	
- 500 slides, packaged as 20 boxes of 25 slides	
Description	POR Cat. Number
White Tab	FT4981IGLPLUS-001

Epredia Superfrost PLUS Gold Slides	
- Unit of Measure - Case (CV)	
- 1,440 slides, packaged as 20 boxes of 72 slides	
Description	POR Cat. Number
White Tab with control box on front	ES-3012-GLPLUS-001



Superfrost PLUS Control Slides

- Designed to run the control and sample on the same slide
- The control sample is denoted by a box at the frosted end of the slide
- Available with PLUS formula to provide the maximum tissue adhesion

Eprelia Superfrost PLUS Control Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
1" White Tab Red control box on back	ES-334PLUS-001
¾" White Tab Red control box on front	ES-329PLUS-001
¾" White Tab Red control box on back	ES-355PLUS-001
1" White Tab Frosted Box on Back	I8446-PLUS-001
¾" White Tab Frosted Box on Back	I8432-PLUS-001




Superfrost Control Slides

- Designed to run the control and sample on the same slide
- The control sample is denoted by a box at the frosted end of the slide

Eprelia Superfrost Control Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
¾" White Tab Red control box on front	ES-325-001
1" White Tab Red control box on back	ES-361-001




Inkjet PLUS Slides

- Designed for printing instruments that use ink and print directly on the tab of the slide
- Available with PLUS formula to provide the maximum tissue adhesion
- Clipped corners to 45°

Eprelia Inkjet PLUS Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number	KZOO Cat. Number
White	4650WH	4650WH




Colormark PLUS Slides

- For use with automated slide etcher systems
- Self-contained packages ready for loading into your slide labeling system
- Simply remove the cover from the slide package and drop the sleeve into the loading bin of the labeling system
- Available with PLUS formula to provide the maximum tissue adhesion

Eprelia Colormark PLUS Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number	KZOO Cat. Number
White/Black	CM-4951WPLUS-001	B9992010
Blue/Dark Blue	CM-5951WPLUS-001	B9992010BL
Pink/Red	CM-7951WPLUS-001	B9992010PK
Yellow/Red	CM-8951WPLUS-001	B9992010YW
Gold/Black	CM-9951WGL-PS-001	B9992010GL
Aqua/Black	CM-9951WA-PLUS-001	B9992010AQ
Orange/Black	CM-9951WBO-PS-001	B9992010BO
Tan/Black	CM-9951WT-PLUS-001	B9992010TN
Lilac/Black	CM-9951WL-PLUS-001	B9992010LV




Colormark Slides

- Designed for use with automated slide etcher systems
- Self-contained packages ready for loading into your slide labeling system
- Simply remove the cover from the slide package and drop the sleeve into the loading bin of the labeling system

Eprelia Colormark Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number	KZOO Cat. Number
White/Black	CM-4951W-001	B9992000
Pink/Red	CM-7951W-001	B9992000PK
Blue/Dark Blue	CM-5951W-001	B9992000BL
Yellow/Red	CM-8951W-001	B9992000YW
White/Black, 25x75x1.2mm	CM-4955W-001	N/A
Aqua/Black	CM-9951WA-001	B9992000AQ
Orange/Black	CM-9951BO-001	B9992000BO
Tan/Black	CM-9951WT-001	B9992000TN
Gold/Black	CM-9951GL-001	B9992000GL
Lilac/Black	CM-9951WL-001	B9992000LV




Plain Slides

- Plain microscope slides are manufactured to the highest quality standards for use in hematology, microbiology and urinalysis as well as other general laboratory and research applications
- All slides are pre-cleaned and feature finely ground edges for easy handling and immediate use
- The surface of each slide is completely plain without a frosted tab

Epredia Plain Slides

- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Plain Slides	2950-001T
Plain Slides	500-001T
Plain Slides 3 x 1"	2950F-001T
w/ clipped corners	2950-4-001



Frosted Slides

- Frosted tab which provides a smooth frosted finish that can be used for manual notation of patient identification or the application of printed labels
- Manufactured using the same high-quality process as the plain slides

Epredia Frosted Slides

- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number	KZOO Cat. Number
Double Frosted, One End	2951D-001T	N/A
Frosted Tab on One Side	4941-8515-001T	N/A
Double Frosted Tab	510D-001T	4801
Fully Frosted on One Side	4958-8516-001T	N/A
Fully Frosted Slides 3 x 1"	2958F-006	N/A



Cytology Slides

- Available in adhesion and non-adhesion (coated and non-coated)
- Adhesion coating applied to the surface of the slide that assists with adhering cytology specimen samples permanently to the surface of the slide

Epredia Ringed Cytoslides™

- Unit of Measure - Case (CV)
- 1,500 slides, packaged as 15 boxes of 100 slides

Description	POR and KZOO Cat. Number
Non-coated circle on back	5991051
Non-coated circle on front	5991059
Coated, circle on back	5991056
DoubleCytoSlide™ Coated	5991055
DoubleCytoSlide	5991054
w/ Black Mask	5991057
w/ Blue Mask	B5991050

Epredia Ringed Cytoslides™


- Unit of Measure - Case (CV)
- 1,350 slides, packaged as 54 boxes of 25 slides

Description	POR Cat. Number	KZOO Cat. Number
Megafunnel Slide, Coated	5991026	N/A

Epredia Ringed Cytoslides™

- Unit of Measure - Carton (CAR)
- 2,000 slides, packaged as 80 boxes of 25 slides

Description	POR and KZOO Cat. Number
Coated, circle on back	B1000731
Double Circle, Coated	5991006



Premium Slides

- Premium plain and frosted slides are manufactured through a unique, proprietary cleaning process intended to provide consistent results and ease of interpretation
- This process allows these slides to be used in hematology applications by delivering consistently crisp, feathered-edge blood smears with good morphology
- The premium frosted clipped-corner slide is specifically designed for use with automated hematology instrumentation.

Epredia Premium Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Superfrost white (no SPECIMEN printed-key)	4941X-006
Superfrost white with four clipped corners	4951W-S-4-4370-006
Superfrost white 3 x 1"	4951FX-006
Plain 3 x 1"	2950FX-006


Epredia Premium Slides
 - Unit of Measure - Case (CV)
 - 3,600 slides, packaged as 50 boxes of 72 slides

Description	POR Cat. Number
Plain	2950WX-003

- Slides are packaged in a moisture-control barrier bag

Epredia Premium Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Plain Slides with humidity control	3011-002A
Two 10mm circles with humidity control	3032-002A
Superfrost White Tab with humidity control	3050-002A
Superfrost White Tab, Epredia Imprint with humidity control	3063-002A



Large Format Slides

- Used for creating whole mount sections
- Available in a variety of sizes and thicknesses

Epredia Large Format Slides
 - Unit of Measure - Case (CV)
 - 720 slides, packaged as 10 boxes of 72 slides

Description	POR Cat. Number
Superfrost PLUS 3" x 2" x 1mm	4957PLUS-006
Superfrost 3 x 2"	4957-006
Superfrost 3 x 2" x 1.2mm	4957FT-006
Plain 3 x 2" x 1.2mm	2957FT-001
Plain 3 x 2" x 1mm	2957WFT-006
Plain 75 x 51 x 1mm	2957-001
Plain 75 x 38 x 1mm	2956-001
Plain 3 x 1.5" x 1mm	2956FT-006

Epredia Large Format Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Plain 25 x 75 x 1.2mm	2954-001
Plain 3 x 1" x 1.2mm	2954F-001
Frosted 25 x 75 x 1.2mm	2955-001
Frosted 3 x 1" x 1.2mm	2955F-001
Frosted 3 x 1" x 1.2mm with SPECIMEN key	4955F-001

Petrographic Slides

- Used to mount thin sections of rocks, minerals, building materials or bone

Epredia Petrographic Slides
 - Unit of Measure - Case (CV)
 - 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Frosted 27x46x1.2mm	2795-001
26x46x1.2mm	2646-001
27x46x1.2mm	2746-001
25x45x1.2mm	2545-001



Specialty Slides

- Localize cells in an easily-identifiable area on the slide surface with EpreDia diagnostic slides

EpreDia Specialty Slides

- Unit of Measure - Case (CV)
- 1,440 slides, packaged as 20 boxes of 72 slides

Description	POR Cat. Number
Frosted Tab with Two 10mm Circles	2960-001
Frosted Tab with Two 15mm Circles	2961-001
Superfrost top justified 13mm Ring Slides with clipped corners	ES-333W-4-72-001
Superfrost bottom justified 13mm Ring Slides with clipped corners	ES-3181W-4-72-001

EpreDia Specialty Slides

- Unit of Measure - Case (CAR)
- 1,500 slides, packaged as 15 boxes of 100 slides

Description	POR Cat. Number
Multi-Spot 3 x 1" with 12 wells (6mm)	9991090
Multi-Spot 3 x 1" with 12 wells (9mm)	9991095


Microarray Slides

- Compatible with standard microarray spotting and detection systems
- Packaged in a secure, no-contact mBox™ and vacuum sealed in a foil barrier bag

EpreDia Microarray Slides


- Unit of Measure - Box (BX)
- 20 slides, packaged as 1 box of 20 slides

Description	POR Cat. Number
SuperChip™ Microarray Slides, Aminopropylsilane, 70 x 20 mm	C18-5131-M20
SuperChip Microarray Slides, Poly-L-Lysine, 70 x 20 mm	C40-5257-M20
UltraClean Microarray Slides, Regular Surface	C22-5128-M20



Coverglass

- Cover slips are manufactured with colorless borosilicate glass. EpreDia cover slips feature high spectral-transmission quality and a finely adaptive refractive index suited for microscopy
- They are available in a variety of sizes, thicknesses and shapes – including squares, rectangles and circles
- All 1oz packages
- All coverglass is sold by the case which includes 10 1oz boxes of coverglass
- Thicknesses range from
 - #0 (0.085 to 0.13 mm)
 - #1 (0.13 to 0.16 mm)
 - #1.5 (0.16 to 0.19 mm)
 - #2 (0.19 to 0.23 mm)



Rectangles

Dimensions	Thickness	POR Cat. Number	KZOO Cat. Number	Approx number of coverslips per package
22 x 30 mm	#1	22X30-1-001G	N/A	119
22 x 30 mm	#1.5	22X30-1.5-001G	N/A	103
22 x 40 mm	#1	22X40-1-001G	102240	89
22 x 40 mm	#1.5	22X40-1.5-001G	N/A	72
22 x 40 mm	#2	22X40-2-001G	N/A	59
22 x 50 mm	#1	22X50-1-001G	102250	71
22 x 50 mm	#1.5	22X50-1.5-001G	152250	58
22 x 60 mm	#1	22X60-1-001G	102260	58
22 x 60 mm	#1.5	22X60-1.5-001G	N/A	48
24 x 30 mm	#1	24X30-1-001G	102430	109
24 x 30 mm	#1.5	24X30-1.5-001G	N/A	94
24 x 40 mm	#1	24X40-1-001G	102440	82
24 x 40 mm	#1.5	24X40-1.5-001G	152440	66
24 x 50 mm	#1	24X50-1-001G	102450	65
24 x 50 mm	#1.5	24X50-1.5-001G	152450	56
24 x 50 mm	#2	24X50-2-001G	202450	43
24 x 55 mm	#1	24X55-1-001G	102455	59
24 x 55 mm	#1.5	N/A	152455	51
24 x 60 mm	#1	24X60-1-001G	102460	54
24 x 60 mm	#1.5	24X60-1.5-001G	152460	47
24 x 60 mm	#2	24X60-2-001G	202460	36

Premium Coverglass


- Premium coverglass from EpreDia is packaged in specialty boxes and includes a desiccant

Dimensions	Thickness	Cat. Number	Approx number of coverslips per package
22 x 40 mm	#1	22X40-1-002G	89
22 x 22 mm	#1	22X22-1-002G	162
24 x 40 mm	#1	24X40-1-002G	82
24 x 50 mm	#1	24X50-1-002G	65
24 x 60 mm	#1	24X60-1-002G	54
48 x 60 mm	#1	48X60-1-002G	27

Specialty Rectangles

- Superslip™ coverslips provide a more slippery surface than traditional coverglass to enhance usability

Dimensions	Thickness	POR Cat. Number	KZOO Cat. Number	Approx number of coverslips per package
24 x 40 mm	#1	24X40-1W-YO-001G	12440S	82
24 x 50 mm	#1	24X50-1W-YO-001G	12450S	54
24 x 50 mm	#1.5	N/A	152450S	48
24 x 55 mm	#1	N/A	12455S	48
24 x 60 mm	#1	24X60-1W-YO-001G	12460S	65




Squares

Dimensions	Thickness	POR Cat. Number	KZOO Cat. Number	Approx number of coverslips per package
18 x 18 mm	#1	18X18-1-001G	N/A	242
18 x 18 mm	#1.5	18X18-1.5-001G	N/A	194
22 x 22 mm	#0	22X22-0-001G	N/A	240
22 x 22 mm	#1	22X22-1-001G	102222	162
22 x 22 mm	#1.5	22X22-1.5-001G	152222	130
22 x 22 mm	#2	22X22-2-001G	N/A	120
24 x 24 mm	#1	24X24-1-001G	102424	136
25 x 25 mm	#1	25X25-1-001G	N/A	121
25 x 25 mm	#2	25X25-2-001G	N/A	92

Circles

Diameter	Thickness	POR Cat. Number	Approx number of coverslips per package
10 mm	#1	10CIR-1-001G	962
12 mm	#1	12CIR-1-001G	685
12 mm	#1.5	12CIR-1.5-001G	567
12 mm	#2	12CIR-2-001G	467
15 mm	#1	15CIR-1-001G	460
15 mm	#1.5	15CIR-1.5-001G	362
15 mm	#2	15CIR-2-001G	307
18 mm	#1	18CIR-1-001G	320
18 mm	#2	18CIR-1.5-001G	213
18 mm	#1.5	18CIR-2-001G	252
22 mm	#1	22CIR-1-001G	199
22 mm	#2	22CIR-1.5-001G	143
22 mm	#1.5	22CIR-2-001G	168
25 mm	#1	25CIR-1-001G	165
25 mm	#1.5	25CIR-1.5-001G	130
25 mm	#2	25CIR-2-001G	110
30 mm	#1.5	30CIR-1.5-001G	91
40 mm	#1	40CIR-1-001G	60
50 mm	#1.5	50CIR-1.5-001G	33



1oz package of coverglass



Package of 10 1oz boxes of coverglass

Coverglass Hoppers

- To be used with EpreDia Clearvue Automated Coverslipper
- 2 hoppers per case
- 575 pieces per hopper for #1 thickness
- 500 pieces per hopper for #1.5 thickness

Dimensions	Thickness	KZOO Cat. Number	Approx number of coverslips per package
24 x 50mm	#1	A79210166	575
24 x 40mm	#1.5	A79210050	500
24 x 50mm	#1.5	A79210051	500
24 x 55mm	#1.5	A79210052	500





Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2025 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. M41092_LTR R042025


A Member of PHC Group