



# Positively Consistent.

Epredia continues to improve glass slides to provide high quality products for the histology lab. Continuous improvement efforts in manufacturing enable Epredia to deliver glass slides with consistent tissue adhesion.

## Confident Staining Starts with Reliable Adhesion

Enhance the quality and consistency of the histology workflow with Epredia's Superfrost PLUS microscope slides. Designed with advanced adhesion properties, these slides help ensure tissue sections remain securely in place during even the most rigorous staining protocols —including H&E, special stains, and immunohistochemistry.

## Take charge of tissue adhesion with Epredia Superfrost PLUS today






## Superfrost PLUS™ Adhesion Slides *+Enhanced+*

**Superfrost PLUS adhesion slides** are designed to secure tissue sections without the need for laboratory applied coatings. Tissue loss is minimised, which results in quick turn-around time and reduced risk.

- Proprietary chemistry of the slide is enhanced to improve adhesion performance stability and reduce impact from environmental conditions.
- Manufactured from extra-white soda-lime glass with very low iron content
- The printed tab ink is impervious to virtually all common laboratory chemicals and reagents.
- Clean, flat, and ready to use
- Manufactured in accordance with US dimensions: 25 x 75 x 1.0 mm
- The EpreDia slide manufacturing facility is ISO13485 certified.
- Can secure frozen and formalin fixed & paraffin-embedded tissue sections.

### Ordering Information

72 pcs per box, 20 boxes per carton.

Product	Cat. No.
<b>EpreDia Superforst Plus <i>+Enhanced+</i></b>	
	Two "+" Markings
White	J1800AMNZ
White with four clipped corners	J7800AMNZ
Blue	J1810AMNZ
Pink	J1820AMNZ
Yellow	J1830AMNZ
Green	J1840AMNZ
Orange	J1850AMNZ
Orange with four clipped corners	J7850AMNZ



Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2025 EpreDia Group. All rights reserved. All trademarks are the property of EpreDia Group and its subsidiaries unless otherwise specified. **M2500063\_A4 R082025**