



Click Top Prefilled Containers

Enhancement with audible click top safety seal

Key enhancements with new prefill container design:

- Single line up between the lid and container to minimize mistreading
- Audible click top safety seal that requires minimal torque
- Dual audible clicks to aide in guaranteeing a fully sealed container
- 95 kPa compliant
- Additional space for sample information
- Item number and lot number now printed on the bottom of the container



Product Description	Cat Number	Qty.	Container Size
---------------------	------------	------	----------------

Epredia Prefilled Containers - Click Top Cutout Cardboard Trays

10% NBF 20mL container	56201	384/cs	20ml
10% NBF 40mL container	56401	216/cs	40ml
10% NBF 60mL container	56601	144/cs	60ml
10% NBF 90mL container	56901	100/cs	90ml
10% NBF 120mL container	561201	75/cs	120ml

Epredia Prefilled Containers - Click Top Bulk Cases

10% NBF 20mL container	59201	192/cs	20ml
10% NBF 40mL container	59401	168/cs	40ml
10% NBF 60mL container	59601	126/cs	60ml
10% NBF 90mL container	59901	120/cs	90ml
10% NBF 120mL container	591201	90/cs	120ml

Epredia Prefilled Containers - Click Top Cardboard Trays

10% NBF 20mL container	LC-0020	784/cs	20ml
10% NBF 40mL container	LC-0040	420/cs	40ml
10% NBF 60mL container	LC-0060	300/cs	60ml
10% NBF 90mL container	LC-0090	240/cs	90ml
10% NBF 120mL container	LC-0120	200/cs	120ml

Epredia Prefilled Containers - Click Top in Handi-Packs

10% NBF 20mL container	51201	192/cs	20ml
10% NBF 40mL container	51401	96/cs	40ml
10% NBF 60mL container	51601	96/cs	60ml
10% NBF 90mL container	51901	96/cs	90ml
10% NBF 120mL container	511201	96/cs	120ml

Epredia Prefilled Containers - Click Top Specimen Transport System (STS) Trays

10% NBF 60mL container	53601	150/cs	60ml
10% NBF 90mL container	53901	125/cs	90ml
10% NBF 120mL container	531201	100/cs	120ml

Visit epredia.com for more information or to view our entire catalogue