



Convenient specimen transport, storage and grossing containers

Epredia Prefilled and Dry Containers



Epredia specimen and grossing containers provide convenient storage and fixation of tissue specimens. The label provides a space for patient information and each prefilled container is marked with a lot number and expiration date. All prefilled containers are labeled with highly visible labels.

Click Top Prefill Containers

Prefilled containers use a click top closure system. They also have a single line up between the lid and container to help minimize the potential for misthreading. The click top safety seal requires minimal torque and provides dual audible clicks to support confidence of a fully sealed container.

Containers are all 95 kPa compliant and made from polystyrene plastic that is completely translucent.



Prefilled Containers with Click Top in Trays

Prefilled containers are separated into cut out cardboard trays. They are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|----------------------------|-------|-------------|-------------------|--------------|
| 56201 Prefilled Container | 20mL | 14mL | 384 | Cutout Trays |
| 56401 Prefilled Container | 40mL | 27mL | 216 | Cutout Trays |
| 56601 Prefilled Container | 60mL | 40mL | 144 | Cutout Trays |
| 56901 Prefilled Container | 90mL | 60mL | 100 | Cutout Trays |
| 561201 Prefilled Container | 120mL | 80mL | 75 | Cutout Trays |



Prefilled Containers with Click Top in Bulk

Prefilled containers packaged in a cardboard box and organized using stackable lids. They are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|----------------------------|-------|-------------|-------------------|-----------|
| 59201 Prefilled Container | 20mL | 14mL | 192 | Bulk |
| 59401 Prefilled Container | 40mL | 27mL | 168 | Bulk |
| 59601 Prefilled Container | 60mL | 40mL | 126 | Bulk |
| 59901 Prefilled Container | 90mL | 60mL | 120 | Bulk |
| 591201 Prefilled Container | 120mL | 80mL | 90 | Bulk |



Prefilled Containers with Click Top in Cardboard Trays

Prefilled containers are packaged in a cardboard trays. They are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|-----------------------------|-------|-------------|-------------------|-----------------|
| LC-0020 Prefilled Container | 20mL | 10mL | 784 | Cardboard Trays |
| LC-0040 Prefilled Container | 40mL | 20mL | 420 | Cardboard Trays |
| LC-0060 Prefilled Container | 60mL | 30mL | 300 | Cardboard Trays |
| LC-0090 Prefilled Container | 90mL | 45mL | 240 | Cardboard Trays |
| LC-0120 Prefilled Container | 120mL | 60mL | 200 | Cardboard Trays |



Prefilled Containers with Click Top in Handi-pack

Prefilled containers are separated into Handi-Packs containers within the package. Prefilled containers are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|----------------------------|-------|-------------|-------------------|------------|
| 51201 Prefilled Container | 20mL | 10mL | 192 | Handi-Pack |
| 51401 Prefilled Container | 40mL | 20mL | 96 | Handi-Pack |
| 51601 Prefilled Container | 60mL | 30mL | 96 | Handi-Pack |
| 51901 Prefilled Container | 90mL | 45mL | 96 | Handi-Pack |
| 511201 Prefilled Container | 120mL | 60mL | 96 | Handi-Pack |



Prefilled Containers with Click Top in Specimen Transport System

Prefilled containers are separated into Specimen Transport System (STS) trays. STS trays have cut out sections for each container to allow for easy transportation. These trays are allowed in some operating suites. Prefilled containers are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|----------------------------|-------|-------------|-------------------|-----------|
| 53601 Prefilled Container | 60mL | 40mL | 150 | STS |
| 53901 Prefilled Container | 90mL | 60mL | 125 | STS |
| 531201 Prefilled Container | 120mL | 80mL | 100 | STS |




Traditional closure prefill containers

The containers are made from polypropylene plastic, which is completely translucent, and utilize a Bio-Tite lid closure system for the cap.

Prefilled Containers in Specimen Transport System

Prefilled containers are separated into Specimen Transport System (STS) trays. STS trays have cut out sections for each container to allow for easy transportation. These trays are allowed in some operating suites. Prefilled containers are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|----------------------------|-------|-------------|-------------------|-----------|
| 53151 Prefilled Container | 15mL | 10mL | 384 | STS |
| 53301 Prefilled Container | 30mL | 20mL | 216 | STS |
| 531801 Prefilled Container | 180mL | 120mL | 64 | STS |
| 534801 Prefilled Container | 480mL | 320mL | 27 | STS |



Prefilled Containers in Bulk

Prefilled containers are packaged in a cardboard box and organized using stackable lids. Prefilled containers are non-sterile and 1/2 or 2/3 filled with 10% Neutral Buffered Formalin (NBF) depending on the configuration.

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|------------------------------|-------|-------------|-------------------|-----------------|
| 591801 Prefilled Container | 180mL | 120mL | 50 | Bulk |
| 594801 Prefilled Container | 480mL | 320mL | 24 | Cardboard Trays |
| 5912001 Prefilled Container | 1.2L | 600mL | 12 | Bulk |
| 599601-9 Prefilled Container | 960mL | 450mL | 9 | Cardboard Trays |



Prefilled Containers in Cardboard Trays

Prefilled containers are packaged in a cardboard trays. They are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF).

| SKU and Description | Size | Fill Volume | Quantity per Case | Packaging |
|--------------------------------|-------|-------------|-------------------|-----------------|
| RBBP-0240 Prefilled Container | 240mL | 120mL | 50 | Cardboard Trays |
| RBBP-0480 Prefilled Container | 480mL | 240mL | 24 | Cardboard Trays |
| RBBP-0960 Prefilled Container | 960mL | 480mL | 9 | Cardboard Trays |
| RBH-0180-S Prefilled Container | 180mL | 90mL | 64 | Cardboard Trays |



Zinc Formalin Prefilled Containers in Bulk

Prefilled containers use a Bio-Tite closure system and are packaged in a cardboard box and organized using stackable lids. They are non-sterile and 2/3 filled with 10% Zinc Buffered Formalin.

| SKU and Description | Size | Fill Volume | Quantity per Case |
|--------------------------------|-------|-------------|-------------------|
| 59201ZF Prefilled Container | 20mL | 14 | 192 |
| 59401ZF Prefilled Container | 40mL | 27 | 168 |
| 59601ZF Prefilled Container | 60mL | 40 | 126 |
| 59901ZF Prefilled Container | 90mL | 60 | 120 |
| 591201ZF Prefilled Container | 120mL | 80 | 90 |
| 591801ZF Prefilled Container | 180mL | 120 | 50 |
| 594801ZF Prefilled Container | 480mL | 320 | 24 |
| 599601ZF-9 Prefilled Container | 960mL | 450 | 9 |



Dry Containers

Epredia dry specimen containers provide convenient materials for fixation and transportation of tissue specimens. Specimen containers are made from polystyrene plastic and are available in clear and opaque options.

Opaque Specimen Containers in Sleeves

The containers are non-sterile and have snap top lids. Containers and lids are packaged in sleeves within the cardboard box.

| SKU and Description | Size | Quantity per Case | Packaging |
|----------------------|-------|-------------------|-----------|
| 5312 Empty Container | 360mL | 350 | Sleeves |
| 5316 Empty Container | 480mL | 350 | Sleeves |
| 5332 Empty Container | 960mL | 250 | Sleeves |
| 5333 Empty Container | 1L | 100 | Sleeves |
| 5364 Empty Container | 2L | 50 | Sleeves |
| 5357 Empty Container | 5.7L | 25 | Sleeves |



Opaque Specimen Containers in Bulk Packs

The containers are non-sterile and have snap top lids.

| SKU and Description | Size | Quantity per Case | Packaging |
|-----------------------|-------|-------------------|-----------|
| 41751 Empty Container | 16oz | 100 | Bulk |
| 41753 Empty Container | 86oz | 25 | Bulk |
| 41754 Empty Container | 172oz | 10 | Bulk |



Clear Specimen Containers in Bulk Packs

Specimen containers are non-sterile and completely translucent.

| SKU and Description | Size | Quantity per Case | Packaging |
|------------------------|-------|-------------------|-----------|
| 53152 Empty Container | 15mL | 384 | Bulk |
| 53302 Empty Container | 30mL | 216 | Bulk |
| 53602 Empty Container | 60mL | 150 | Bulk |
| 53902 Empty Container | 90mL | 125 | Bulk |
| 531202 Empty Container | 120mL | 100 | Bulk |
| 531802 Empty Container | 180mL | 150 | Bulk |
| 534802 Empty Container | 480mL | 60 | Bulk |



Clear Specimen Containers in Bulk Packs with Stackable Lids

The containers are non-sterile, completely translucent and have lids that are conducive for stacking multiple containers on top of one another.

| SKU and Description | Size | Quantity per Case | Packaging |
|-------------------------|-------|-------------------|-----------|
| 59202 Empty Container | 20mL | 384 | Bulk |
| 59402 Empty Container | 40mL | 216 | Bulk |
| 59602 Empty Container | 60mL | 150 | Bulk |
| 59902 Empty Container | 90mL | 125 | Bulk |
| 591202 Empty Container | 120mL | 100 | Bulk |
| 599602 Empty Container | 960mL | 12 | Bulk |
| 5912002 Empty Container | 1.2L | 12 | Bulk |



Specialty Prefilled Containers

Prefill containers containing a variety of different solutions. Applications for the alternative solutions can vary from cytology applications or specific histology needs.

- **Formal-Fixx**

- Phosphate-free
- Use for storage of all tissue types
- Stores specimens without excessive shrinking or hardening
- Available premixed, ready-to-use



| SKU and Description | Contents | Size | Fill Volume | Quantity per Case |
|-----------------------------|-------------|-------|-------------|-------------------|
| 9990914 Prefilled Container | Formal-Fixx | 60mL | 30 | 125 |
| 9990916 Prefilled Container | Formal-Fixx | 90mL | 45 | 125 |
| 9990918 Prefilled Container | Formal-Fixx | 120mL | 60 | 125 |



Prostate Kit with Prefilled Containers

In addition to the 16 vials and biohazard bag, the Prostate Kit includes instructions printed on the inside of the box to help prepare the prostate sample for transport.

| SKU and Description | Contents | Size | Quantity per Case |
|----------------------------|----------|------|-------------------|
| 29016 Prostate Kit Prefill | 10% NBF | 20mL | 16 |



To find out about all **laboratory solutions** from Epredia, contact your **Epredia Account Executive**.

To place an **order**, contact **Epredia Customer Care** at ep-cs-kalamazoo@epredia.com; (800) 522-7270.

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2025 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. **M2400058_US R102025**

