



Convenient specimen transport, storage and grossing containers

Epredia Prefilled and Dry Containers



Epredia specimen and grossing containers provide convenient storage and fixation of tissue specimens. The label provides a space for patient information and each prefilled container is marked with a lot number and expiration date. All prefilled containers are labeled with highly visible labels.

Click Top Prefill Containers

Prefilled containers use a click top closure system. They also have a single line up between the lid and container to help minimize the potential mistreading. The click top safety seal requires minimal torque and provides dual audible clicks to support confidence of a fully sealed container.

Containers are all 95 kPa compliant and made from polystyrene plastic that is completely translucent.



Prefilled Containers with Click Top in Trays

Prefilled containers are separated into cut out cardboard trays. They are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
56201 Prefilled Container	20mL	14mL	384	Cutout Trays
56401 Prefilled Container	40mL	27mL	216	Cutout Trays
56601 Prefilled Container	60mL	40mL	144	Cutout Trays
56901 Prefilled Container	90mL	60mL	100	Cutout Trays
561201 Prefilled Container	120mL	80mL	75	Cutout Trays



Prefilled Containers with Click Top in Bulk

Prefilled containers packaged in a cardboard box and organized using stackable lids. They are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
59201 Prefilled Container	20mL	14mL	192	Bulk
59401 Prefilled Container	40mL	27mL	168	Bulk
59601 Prefilled Container	60mL	40mL	126	Bulk
59901 Prefilled Container	90mL	60mL	120	Bulk
591201 Prefilled Container	120mL	80mL	90	Bulk



Prefilled Containers with Click Top in Cardboard Trays

Prefilled containers are packaged in a cardboard trays. They are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF).

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
LC-0020 Prefilled Container	20mL	10mL	784	Cardboard Trays
LC-0040 Prefilled Container	40mL	20mL	420	Cardboard Trays
LC-0060 Prefilled Container	60mL	30mL	300	Cardboard Trays
LC-0090 Prefilled Container	90mL	45mL	240	Cardboard Trays
LC-0120 Prefilled Container	120mL	60mL	200	Cardboard Trays



Traditional closure prefill containers:

The containers are made from polypropylene plastic, which is completely translucent, and utilize a Bio-Tite lid closure system for the cap.



Prefilled Containers in Handi-pack

Prefilled containers are separated into Handi-Packs containers within the package. Prefilled containers are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF). All prefilled containers are labeled with black and white labels.


SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
51201 Prefilled Container	20mL	10mL	192	Handi-Pack
51401 Prefilled Container	40mL	20mL	96	Handi-Pack
51601 Prefilled Container	60mL	30mL	96	Handi-Pack
51901 Prefilled Container	90mL	45mL	96	Handi-Pack
511201 Prefilled Container	120mL	60mL	96	Handi-Pack



Prefilled Containers in Specimen Transport System

Prefilled containers are separated into Specimen Transport System (STS) trays. STS trays have cut out sections for each container to allow for easy transportation. These trays are allowed in some operating suites. Prefilled containers are non-sterile and 2/3 filled with 10% Neutral Buffered Formalin (NBF).

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
53151 Prefilled Container	15mL	10mL	384	STS
53301 Prefilled Container	30mL	20mL	216	STS
53601 Prefilled Container	60mL	40mL	150	STS
53901 Prefilled Container	90mL	60mL	125	STS
531201 Prefilled Container	120mL	80mL	100	STS
531801 Prefilled Container	180mL	120mL	64	STS
534801 Prefilled Container	480mL	320mL	27	STS



Prefilled Containers in Bulk

Prefilled containers are packaged in a cardboard box and organized using stackable lids. Prefilled containers are non-sterile and 1/2 or 2/3 filled with 10% Neutral Buffered Formalin (NBF) depending on the configuration.

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
591801 Prefilled Container	180mL	120mL	50	Bulk
594801 Prefilled Container	480mL	320mL	24	Cardboard Trays
5912001 Prefilled Container	1.2L	600mL	12	Bulk
599601-9 Prefilled Container	960mL	450mL	9	Cardboard Trays



Prefilled Containers in Cardboard Trays

Prefilled containers are packaged in a cardboard trays. They are non-sterile and 1/2 filled with 10% Neutral Buffered Formalin (NBF).

SKU and Description	Size	Fill Volume	Quantity per Case	Packaging
RBBP-0240 Prefilled Container	240mL	120mL	50	Cardboard Trays
RBBP-0480 Prefilled Container	480mL	240mL	24	Cardboard Trays
RBBP-0960 Prefilled Container	960mL	480mL	9	Cardboard Trays
RBH-0180-S Prefilled Container	180mL	90mL	64	Cardboard Trays



Zinc Formalin Prefilled Containers in Bulk

Prefilled containers use a Bio-Tite closure system and are packaged in a cardboard box and organized using stackable lids. They are non-sterile and 2/3 filled with 10% Zinc Buffered Formalin.

SKU and Description	Size	Fill Volume	Quantity per Case
59201ZF Prefilled Container	20mL	14	192
59401ZF Prefilled Container	40mL	27	168
59601ZF Prefilled Container	60mL	40	126
59901ZF Prefilled Container	90mL	60	120
591201ZF Prefilled Container	120mL	80	90
591801ZF Prefilled Container	180mL	120	50
594801ZF Prefilled Container	480mL	320	24
599601ZF-9 Prefilled Container	960mL	450	9



Dry Containers

Epreidia dry specimen containers provide convenient materials for fixation and transportation of tissue specimens. Specimen containers are made from polystyrene plastic and are available in clear and opaque options.

Opaque Specimen Containers in Sleeves

The containers are non-sterile and have snap top lids. Containers and lids are packaged in sleeves within the cardboard box.

SKU and Description	Size	Quantity per Case	Packaging
5312 Empty Container	360mL	350	Sleeves
5316 Empty Container	480mL	350	Sleeves
5332 Empty Container	960mL	250	Sleeves
5333 Empty Container	1L	100	Sleeves
5364 Empty Container	2L	50	Sleeves
5357 Empty Container	5.7L	25	Sleeves



Opaque Specimen Containers in Bulk Packs

The containers are non-sterile and have snap top lids.

SKU and Description	Size	Quantity per Case	Packaging
41751 Empty Container	16oz	100	Bulk
41753 Empty Container	86oz	25	Bulk
41754 Empty Container	172oz	10	Bulk



Clear Specimen Containers in Bulk Packs

Specimen containers are non-sterile and completely translucent.

SKU and Description	Size	Quantity per Case	Packaging
53152 Empty Container	15mL	384	Bulk
53302 Empty Container	30mL	216	Bulk
53602 Empty Container	60mL	150	Bulk
53902 Empty Container	90mL	125	Bulk
531202 Empty Container	120mL	100	Bulk
531802 Empty Container	180mL	150	Bulk
534802 Empty Container	480mL	60	Bulk



Clear Specimen Containers in Bulk Packs with Stackable Lids

The containers are non-sterile, completely translucent and have lids that are conducive for stacking multiple containers on top of one another.

SKU and Description	Size	Quantity per Case	Packaging
59202 Empty Container	20mL	384	Bulk
59402 Empty Container	40mL	216	Bulk
59602 Empty Container	60mL	150	Bulk
59902 Empty Container	90mL	125	Bulk
591202 Empty Container	120mL	100	Bulk
599602 Empty Container	960mL	12	Bulk
5912002 Empty Container	1.2L	12	Bulk



Specialty Prefilled Containers

Prefill containers containing a variety of different solutions. Applications for the alternative solutions can vary from cytology applications or specific histology needs.

- **Formal-Fixx**

- Phosphate-free
- Use for storage of all tissue types
- Stores specimens without excessive shrinking or hardening
- Available premixed, ready-to-use

- **Glyo-Fixx**

- Noncarcinogenic, nonsensitizing, and environmentally friendly substitute for 10% NBF
- Produces distinct and crisp cellular morphology
- Fixes most antigens without impairing immunoreactivity
- Allows special stain dyes, including silver procedures, to be sharp without background staining



SKU and Description	Contents	Size	Fill Volume	Quantity per Case
9990914 Prefilled Container	Formal-Fixx	60mL	30	125
9990916 Prefilled Container	Formal-Fixx	90mL	45	125
9990918 Prefilled Container	Formal-Fixx	120mL	60	125
9990920 Prefilled Container	Glyo-Fixx	20mL	10	180
9990924 Prefilled Container	Glyo-Fixx	60mL	30	125
9990928 Prefilled Container	Glyo-Fixx	120mL	125	60



Prostate Kit with Prefilled Containers

In addition to the 16 vials and biohazard bag, the Prostate Kit includes instructions printed on the inside of the box to help prepare the prostate sample for transport.

SKU and Description	Contents	Size	Quantity per Case
29016 Prostate Kit Prefill	10% NBF	20mL	16



To find out about all laboratory solutions from Epredia, contact your **Epredia Account Executive** or **Epredia Customer Care**.

Customer Care

ep-cs-kalamazoo@epredia.com

(800) 522-7270

Find out more at [epredia.com](https://www.epredia.com)