



# When a touch of pink helps to see clearer

**EpreDia introduces PINK-Fix for easy visualization of sample.**

Once tissue samples are removed from the body, starting the fixation process immediately is necessary, specifically when working with small biopsies. Rose colored formalin in pre-filled containers is designed to visualize small, tissue fragments and biopsy-sized specimens better.

- The colored solution improves the visibility of grey or opaque tissue fragments, aiding grossing technicians and pathologists throughout the grossing process.
- Allowing technicians to visualize small biopsies at the bench may provide better efficiency and mitigates the error of tissue loss during processing and embedding.
- Further visualization of small tissue biopsies throughout the workflow process helps to further maintain the integrity of the tissue.

When looking at small tissue biopsies, PINK-Fix will help you make a difference in your workflow!

Item Number	Description	Order Size	Certification
5735P-LP	Fixative, 10% PINK-Fix	4 x 5 L	IVD
51201P-LP-Y	Prefill, 10% PINK-Fix 20/10 ml, yellow cap	24 x 10 ml	IVD
51601P-LP-Y	Prefill, 10% PINK-Fix 60/30 ml, yellow cap	35 x 30 ml	IVD
531251P-LP-Y	Prefill, 10% PINK-Fix 125/60 ml, yellow cap	24 x 60 ml	IVD

Please contact your **EpreDia Sales Representative** for further information.



Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2024 EpreDia Group. All rights reserved. All trademarks are the property of EpreDia Group and its subsidiaries unless otherwise specified. **M2400098 R092024**

  
A Member of PHC Group