



Enhance traceability

Comprehensive sample identification demands advanced printing capabilities.

epredia™ VEGA laser cassette printers offer:

- High-resolution, durable printing
- Up to 100-character 2D barcodes
- Fast, five seconds per cassette printing
- Intelligent cassette selection from six-hopper carousel
- Quiet printing, < 65 dB



**epredia**



- High-resolution, 600-dpi barcode to include more information
- High-quality barcode printing may reduce instances of reprints and non-scannable barcodes



- LIS connectivity to streamline workflow
- Easy-to-use connection for smooth installation



- Reduce operational cost with no printhead replacements
- Robust laser technology consistently prints within five seconds for a faster turnaround time



Ordering Information

	Epredia Ref.
Epredia VEGA Laser Cassette Printer	
Epredia VEGA laser cassette printer with all-in-one PC, six hoppers and digital imager	VEGA05-FT-PCDI
Epredia VEGA laser cassette printer with 6 hoppers	8000038
Epredia VEGA laser cassette printer with 6 hoppers with sorter	8000040
Epredia VEGA laser cassette printer with 6 hoppers with DUO	8000041
Epredia VEGA laser cassette printer with 8 hoppers	8000062
Accessories	
Side delivery system for VEGA	9001140
Cassette sorter V2 for VEGA	8000042
Cassette collector for VEGA sorter, 4 x 25 cassette capacity	9001387
Duo unit V2 for VEGA	8000043
Trays couple for VEGA Duo sorter, 2 x 10 cassette capacity	9001547
Hopper carousel with 6 hoppers, 480 cassette capacity	9001327
Hopper carousel with 8 hoppers, 640 cassette capacity	9001521
13.3 in. All-in-one PC	B8131005MD



For more information, contact your local Epredia Account Executive or visit [epredia.com](https://www.epredia.com)

