



SlideMate™ Laser

Enter a new era of precision printing using laser technology

Clarity

- High resolution, 600dpi printing
- Laser-optimized slides provide clear and robust printing surface
- Adhesion slides use same technology as Superfrost PLUS™ slides

Consistency

- True on-demand slide barcode printing
- Reliable laser technology may reduce instances of instrument downtime
- No additional consumables required such as ink, bulbs or ribbon

Confidence

- Improved barcode quality can enhance downstream workflow with improved scannability
- Installation and connectivity support provided by experienced Epreidia team

Epreidia Ref.	
Instruments	
SlideMate Laser w/ Filters	8000048A
SlideMate Laser + Total Care Package and Filters	8000048B

Accessories	Quantity	Epreidia Ref.
SlideMate Laser Hopper	2	9001644
SlideMate Laser Chute	2	9001645
SlideMate Laser filters	5	4000850
SlideMate Laser brush	1	20000001
SlideMate Laser stylus	1	20000002





SlideMate Laser Slides

Achieve optimal printing performance

Epredia's SlideMate Laser slides are specifically designed for compatibility with SlideMate Laser.

Key features of these laser optimized slides include:

- 25x75x1mm standard slide size
- Available in adhesion/non adhesion formulations
- Adhesion technology identical to Superfrost PLUS
- Consistent optical clarity based on established Epredia manufacturing procedures
- Available in 5 colors, white, blue, pink, yellow, and green

Key benefits from SlideMate Laser slides include:

- Indelible high resolution printing
- Consistent legible labelling on slide tab
- High success rate of barcode scanning on downstream instrumentation
- Confidently adhere tissue samples from microtomy to sample analysis

SlideMate Laser PLUS Slides	Epredia Ref.
White	LS-40411-PS8523-1
Blue	LS-50411-PS8523-1
Green	LS-60411-PS8523-1
Pink	LS-70411-PS8523-1
Yellow	LS-80411-PS8523-1

SlideMate Laser Non-adhesion Slides	Epredia Ref.
White	LS-40411-8524-1
Blue	LS-50411-8524-1
Green	LS-60411-8524-1
Pink	LS-70411-8524-1
Yellow	LS-80411-8524-1

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2024 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. **M41093_LTR R022024**

