



SlideMate™ Laser

Enter a new era of precision printing using laser technology

Clarity

- High resolution, 600dpi printing
- Laser-optimized slides provide clear and robust printing surface
- Adhesion slides use same technology as Superfrost PLUS™ slides

Consistency

- True on-demand slide barcode printing
- Reliable laser technology may reduce instances of instrument downtime
- No additional consumables required such as ink, bulbs or ribbon

Confidence

- Improved barcode quality can enhance downstream workflow with improved scannability
- Installation and connectivity support provided by experienced EpreDia team

		EpreDia Ref.
Instruments		
SlideMate Laser Printer		8000048
Accessories		
	Quantity	
SlideMate Laser Hopper	2	9001644
SlideMate Laser Chute	2	9001645
SlideMate Laser filters	5	4000850
SlideMate Laser brush	1	20000001
SlideMate Laser stylus	1	20000002





SlideMate Laser Slides

Achieve optimal printing performance

Epredia's SlideMate Laser slides are specifically designed for compatibility with SlideMate Laser.

Key features of these laser optimized slides include:

- 25x75x1mm standard slide size
- Available in adhesion/non adhesion formulations
- Adhesion technology identical to Superfrost PLUS
- Consistent optical clarity based on established Epredia manufacturing procedures
- 45° angle safety corners slide available for downstream instrumentation

Key benefits from SlideMate Laser slides include:

- Indelible high resolution printing
- Consistent legible labelling on slide tab
- Achieved 100% barcode scannability with the Epredia E1000 Dx Digital Pathology Solution & the 3DHISTECH PANNORAMIC 1000 DX Digital Pathology Scanner
- Confidently adhere tissue samples from microtomy to sample analysis

SlideMate Laser PLUS Slides	Epredia Ref.
White	LS-4041IPS8523-1CE
Blue	LS-5041IPS8523-1CE
Green	LS-6041IPS8523-1CE
Pink	LS-7041IPS8523-1CE
Yellow	LS-8041IPS8523-1CE
White, safety corners (clipped)	LS4041PS8566-4-1CE
Blue, safety corners (clipped)	LS5041PS8566-4-1CE
Green, safety corners (clipped)	LS6041PS8566-4-1CE
Pink, safety corners (clipped)	LS7041PS8566-4-1CE
Yellow, safety corners (clipped)	LS8041PS8566-4-1CE

SlideMate Laser Slides (non-adhesion)	Epredia Ref.
White	LS-4041I-8524-1CE
Blue	LS-5041I-8524-1CE
Green	LS-6041I-8524-1CE
Pink	LS-7041I-8524-1CE
Yellow	LS-8041I-8524-1CE
White, safety corners (clipped)	LS-4041I-4-8560C
Blue, safety corners (clipped)	LS-5041I-4-8560C
Green, safety corners (clipped)	LS-6041I-4-8560C
Pink, safety corners (clipped)	LS-7041I-4-8560C
Yellow, safety corners (clipped)	LS-8041I-4-8560C

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2025 Epredia Group. All rights reserved. All trademarks are the property of Epredia and its subsidiaries unless otherwise specified. Please reach out to your local sales representative for more information. **M2500008_EMEA_A4_R122025**

