



SlideMate™ Laser

Enter a new era of precision printing using laser technology

Clarity

- High resolution, 600dpi printing
- Laser-optimized slides provide clear and robust printing surface
- Adhesion slides use same technology as Superfrost PLUS™ slides

Consistency

- True on-demand slide barcode printing
- Reliable laser technology may reduce instances of instrument downtime
- No additional consumables required such as ink, bulbs or ribbon

Confidence

- Improved barcode quality can enhance downstream workflow with improved scannability
- Installation and connectivity support provided by experienced Epredia team

		Epredia Ref.
Instruments		
SlideMate Laser Printer		8000048
Accessories		
	Quantity	
SlideMate Laser Hopper	2	9001644
SlideMate Laser Chute	2	9001645
SlideMate Laser filters	5	4000850
SlideMate Laser brush	1	20000001
SlideMate Laser stylus	1	20000002



SlideMate Laser Slides

Achieve optimal printing performance

Epredia's SlideMate Laser slides are specifically designed for compatibility with SlideMate Laser.

Key features of these laser optimized slides include:

- 25x75x1mm standard slide size
- Available in adhesion/non adhesion formulations
- Adhesion technology identical to Superfrost PLUS
- Consistent optical clarity based on established Epredia manufacturing procedures
- Available in 5 colors, white, blue, pink, yellow, and green

Key benefits from SlideMate Laser slides include:

- Indelible high resolution printing
- Consistent legible labelling on slide tab
- High success rate of barcode scanning on downstream instrumentation
- Confidently adhere tissue samples from microtomy to sample analysis

SlideMate Laser PLUS Slides	Epredia Ref.
White	LS-4041IPS8523-1CE
Blue	LS-5041IPS8523-1CE
Green	LS-6041IPS8523-1CE
Pink	LS-7041IPS8523-1CE
Yellow	LS-8041IPS8523-1CE

SlideMate Laser Non-adhesion Slides	Epredia Ref.
White	LS-4041I-8524-1CE
Blue	LS-5041I-8524-1CE
Green	LS-6041I-8524-1CE
Pink	LS-7041I-8524-1CE
Yellow	LS-8041I-8524-1CE

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2023 Epredia Group. All rights reserved. All trademarks are the property of Epredia and its subsidiaries unless otherwise specified. Please reach out to your local sales representative for more information. M41093_A4_EMEA R122023



POWERED BY SHANDON MICROM MENZEL-GLASER