



SlideMate TSP

Modernized thermal transfer slide printing solution
compatible with a wide range of slide types



Protect specimens and slides

High quality, direct to slide printing

- Print directly on slides and help minimize sample traceability errors due to label transfer
- Lateral slide movement to help reduce slide chipping
- Print on control slides without contacting tissue sample

Enclosed slide hoppers help protect slides from contamination

- Help reduce contamination from laboratory dust and debris by using enclosed hopper solution

Advanced service plans protect your investment

- Comprehensive service plans available tailored to protect your slide printer assets.



Practical slide printing for the modern laboratory

Compatible with a wide range of slide types

- Labs can validate a variety of slide types
- Thermal transfer technology allows printing on large selection of slide types

LCD variable pitch

- Screen pitch adjustment for enhanced ergonomics
- Initialization in <75% less time than other EpreDia thermal transfer slide printers*

Small footprint enables placement at the microtome

- Transition to on-demand printing directly during microtomy


*Data on file

Proven expertise in Slide Printing technology and support

Highly reliable microscope slides

- Validated with SlideMate, SlideMate PLUS, SlideMate Laser and SlideMate Laser PLUS slides
- PLUS slides use the same proven adhesion technology as Superfrost PLUS

Slide Printing and manufacturing expertise worldwide



Reliable service and support from industry experts


- Manufacturer trained service teams



1 8" touch screen



2 Slide Input/Output Channel




3 Barcode reader



4 Slide collection chute (13 max)



5 LED indicators



6 Validated consumables

Instrument Product Ordering Information

Description	Quantity	Cat. Number
SlideMate TSP	1 each	8000066
Hopper	2 per pack	9001644
Chute	1 each	9001645
Ribbon	2 per pack	22000001
Cleaning Kit/Swabs	4 per pack	22000002

Validated Microscope Slide Product Ordering Information

All slides available in quantities of 1 case, 20 boxes of 72 slides, total of 1,440 slides per unit

Description	Cat. Number
SlideMate Laser PLUS White	LS-4041IPS8523-1CE
SlideMate Laser PLUS Blue	LS-5041IPS8523-1CE
SlideMate Laser PLUS Green	LS-6041IPS8523-1CE
SlideMate Laser PLUS Pink	LS-7041IPS8523-1CE
SlideMate Laser PLUS Yellow	LS-8041IPS8523-1CE
SlideMate Laser White	LS-4041I-8524-1CE
SlideMate Laser Blue	LS-5041I-8524-1CE
SlideMate Laser Green	LS-6041I-8524-1CE
SlideMate Laser Pink	LS-7041I-8524-1CE
SlideMate Laser Yellow	LS-8041I-8524-1CE
SlideMate Laser PLUS White w/ Safety Corners	LS4041PS8566-4-1CE
SlideMate Laser PLUS Blue w/ Safety Corners	LS5041PS8566-4-1CE
SlideMate Laser PLUS Green w/ Safety Corners	LS6041PS8566-4-1CE
SlideMate Laser PLUS Pink w/ Safety Corners	LS7041PS8566-4-1CE
SlideMate Laser PLUS Yellow w/ Safety Corners	LS8041PS8566-4-1CE
SlideMate Laser White w/ Safety Corners	LS-4041I-4-8560C
SlideMate Laser Blue w/ Safety Corners	LS-4041I-4-8560C
SlideMate Laser Green w/ Safety Corners	LS-5041I-4-8560C
SlideMate Laser Pink w/ Safety Corners	LS-6041I-4-8560C
SlideMate Laser Yellow w/ Safety Corners	LS-7041I-4-8560C

Description	Cat. Number
SlideMate PLUS White	TT-40418218-PS-W
SlideMate PLUS Blue	TT-50418218-PS-B
SlideMate PLUS Green	TT-60418218-PS-G
SlideMate PLUS Pink	TT-70418218-PS-P
SlideMate PLUS Yellow	TT-80418218-PS-Y
SlideMate White	TT-4041-8217-W
SlideMate Blue	TT-50418217-B
SlideMate Green	TT-60418217-G
SlideMate Pink	TT-70418217-P
SlideMate Yellow	TT-80418217-Y

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2026 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. **M2600008_EMEA_R042026**

