



SlideMate TSP

Modernized thermal transfer slide printing solution compatible with a wide range of slide types



Protect specimens and slides

High quality, direct to slide printing

- Print directly on slides and help minimize sample traceability errors due to label transfer
- Lateral slide movement to help reduce slide chipping
- Print on control slides without contacting tissue sample

Enclosed slide hoppers help protect slides from contamination

- Help reduce contamination from laboratory dust and debris by using enclosed hopper solution

Advanced service plans protect your investment


- Total Care Package and SWAP service plans available to protect your slide printer assets

Proven expertise in Slide Printing technology and support

Highly reliable microscope slides

- Validated with SlideMate, SlideMate PLUS, SlideMate Laser and SlideMate Laser PLUS slides
- PLUS slides use the same proven adhesion technology as Superfrost PLUS

Slide Printing and manufacturing expertise worldwide



Reliable service and support from industry experts

- Manufacturer trained service teams
- SWAP printer exchanged within 24 hours (if reported before 3pm)



Practical slide printing for the modern laboratory

Compatible with a wide range of slide types

- Labs can validate a variety of slide types
- Thermal transfer technology allows printing on large selection of slide types

LCD variable pitch

- Screen pitch adjustment for enhanced ergonomics
- Initialization in <75% less time than other Epredia thermal transfer slide printers*


Small footprint enables placement at the microtome

- Transition to on-demand printing directly during microtomy


*Data on file




1 8" touch screen



2 Slide Input/Output Channel




3 Barcode reader



4 Slide collection chute (13 max)



5 LED indicators



6 Validated consumables



6

Instrument Product Ordering Information

Description	Quantity	Cat. Number
SlideMate TSP	1 each	8000066
Hopper	2 per pack	9001644
Chute	1 each	9001645
Ribbon	2 per pack	22000001
Cleaning Kit/Swabs	4 per pack	22000002
SlideMate TSP with Total Care Package	1 year of warranty, 4 years of SWAP service	8000066B
SWAP Service Plan	1 year of SWAP service	SC-SMTSW0001
Windows 11 Tablet and Docking Station	1 each	A86130057
Keyboard	1 each	AP18268
Mouse	1 each	B81310019
Barcode Scanner w/ Stand	1 each	B35310050

Validated Microscope Slide Product Ordering Information

All slides available in quantities of 1 case, 20 boxes of 72 slides, total of 1,440 slides per unit

Description	Cat. Number
SlideMate PLUS White	TT-4041I-PS-8218-1
SlideMate PLUS Blue	TT-5041I-PS-8218-1
SlideMate PLUS Green	TT-6041I-PS-8218-1
SlideMate PLUS Pink	TT-7041I-PS-8218-1
SlideMate PLUS Yellow	TT-8041I-PS-8218-1
SlideMate White	TT-4041I-8217-001
SlideMate Blue	TT-5041I-8217-001
SlideMate Green	TT-6041I-8217-001
SlideMate Pink	TT-7041I-8217-001
SlideMate Yellow	TT-8041I-8217-001

Description	Cat. Number
SlideMate Laser PLUS White	LS-4041I-PS8523-1
SlideMate Laser PLUS Blue	LS-5041I-PS8523-1
SlideMate Laser PLUS Green	LS-6041I-PS8523-1
SlideMate Laser PLUS Pink	LS-7041I-PS8523-1
SlideMate Laser PLUS Yellow	LS-8041I-PS8523-1
SlideMate Laser White	LS-4041I-8524-1
SlideMate Laser Blue	LS-5041I-8524-1
SlideMate Laser Green	LS-6041I-8524-1
SlideMate Laser Pink	LS-7041I-8524-1
SlideMate Laser Yellow	LS-8041I-8524-1

Find out more at [epredia.com](https://www.epredia.com)

4481 Campus Drive | Kalamazoo, MI 49008 | United States | +1 (800) 522-7270

© 2026 Epredia Group. All rights reserved. All trademarks are the property of Epredia Group and its subsidiaries unless otherwise specified. **M2600008 R042026**

