



Syntri Safeguard Consumables

Diagnostic confidence in your tissue processing hinges upon standardization and consistency throughout the laboratory.

The **Revos rotational tissue processor** tracks and logs reagent usage when paired with reagents featuring Syntri Safeguard RFID technology for the total tissue processing solution.

A Smarter Solution with **Syntri safeguard consumables** allowing the Revos to track every bottle loaded, protocol executed, lot number, and expiration date – helping to reduce downtime and facilitate your workload.





For more information about the Revos tissue processor, please visit epredia.com/revos

Syntri Safeguard Consumables

Description	Cat. Number
Syntri Safeguard™ 10% Neutral Buffered Formalin, 4 x 5 L	5735SS
Syntri Safeguard™ Pen-Fix™, 4 x 5 L	6115SS
Syntri Safeguard Dehydrant 100%, 4 x 5 L	6215SS
Syntri Safeguard Xylene, 4 x 5 L	6615SS
Syntri Safeguard™ Clear-Rite™ 3, 4 x 5 L	6915SS
Syntri Safeguard™ Flex™ 100, 4 x 5 L	8115SS
Syntri Safeguard Isopropyl Alcohol, 4 x 5 L	9517SS
Syntri Safeguard™ Formal-Fixx™, 4 x 5 L	6764215SS
Shandon™ Syntri Safeguard™ Xylene Substitute, 4 x 5 L	6764515SS

 Epredia
4481 Campus Drive
Kalamazoo, MI 49008
UNITED STATES
+1 (800) 522-7270

Find out more at epredia.com

 Enhancing precision
cancer diagnostics