



# Tissue processor preventive maintenance

Today's laboratories are measured by the efficiency and accuracy of their workflow processes. If even a single instrument goes down or malfunctions, patient test results may be at risk.

An investment in an EpreDia Preventative Maintenance Contract is simple, cost-effective way to reduce the chance of your laboratory facing the pitfalls of downtime.

Our Field Service Team will execute all recommended service at the proper intervals. They will thoroughly inspect your tissue processing instrument and perform any needed maintenance or repairs using only genuine OEM-approved parts.

Be confident that your instruments will be up and running when you need them. Contact your local EpreDia representative and learn more about our service offerings today.

Service included	
Electrical and mechanical safety tests	✓
Update software	✓
Investigate and resolve errors identified in the engineer's log	✓
Fluid Transfer System. Adjust and calibrate the following: flush reagents, max fluid / base temperatures, specific gravity meter	✓
Calibrate screen	✓
Verify / update current modification strike status	✓
Inspect the reaction chamber serial number and issue number in production services and update if required	✓
Inspect / repair lid seal	✓
Verify down-draft fan performance	✓
Inspect and clean level sensors	✓
Replace ports, seal, screws and bearings	✓
Test battery back-up system, replace batteries per maintenance schedule	✓
Clean and verify vacuum and pump performance	✓
Replace reagent tubes and fittings as necessary	✓
Inspect, clean, and / or replace on-board reagent bottles	✓
Clean filter baffle and change charcoal and permanganate filters if necessary	✓
Inspect all currents and voltages, adjust / repair	✓
Verify all switches and safety devices are operational	✓
Verify instrument performance	✓

Find out more at [epredia.com](http://epredia.com)

4481 Campus Drive  
Kalamazoo, MI 49008  
UNITED STATES  
+1 (800) 522-7270

Tudor Road, Manor Park  
Runcorn, WA7 1TA  
UNITED KINGDOM  
+44 (0) 800 018 9396  
+44 (0) 1928 534 000

**epredia**  Enhancing precision  
cancer diagnostics

POWERED BY SHANDON  MICRON  MENZEL-GLASER 